

Anti-Anopheles gambiae (African malaria mosquito) AgCactus Antibody Cy3 Conjugated

Catalog Number: DZ41846-Cy3

Description

| | |
|---------------------|--------------------------------------------------------------------------------------------------------------|
| Size | 200 ul/vial |
| Form | Liquid |
| Purification | Immunogen affinity purified. |
| Contents | Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% NaN ₃ . |
| Conjugate | Cy3 Excitation Wavelength: 554 nm Emission Wavelength: 568 nm |
| Suggested Dilutions | Flow Cytometry, Optimal dilutions should be determined by end users. |

Storage

At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Anopheles gambiae (African malaria mosquito) AgCactus Antibody - Cy3

For Research Use Only. Not for use in diagnostic procedures.