

Anti-Gallus gallus (Chicken) SERPINE3 Antibody

Catalog Number: DZ41868-carrier-free

Overview

Product Name	Anti-Gallus gallus (Chicken) SERPINE3 Antibody
Reactive Species	Chicken
Description	Boster Bio Anti-Gallus gallus (Chicken) SERPINE3 Antibody catalog # DZ41868. This antibody reacts with Chicken.
Application	IHC
Clonality	Polyclonal
Formulation	Each vial contains 20mM PBS, 50% glycerol, 0.02% NaN3.
Storage Instructions	At -20°C for one year, at 4°C for one month. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	A0A8V0XIH6

Technical Details

Form	Liquid
Concentration	500 ug/ml
Purification	Immunogen affinity purified.
Suggested Dilutions	Immunohistochemistry (Paraffin-embedded Section), 2-5 ug/ml

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Gallus gallus (Chicken) SERPINE3 Antibody

For Research Use Only. Not for use in diagnostic procedures.