

## Anti-Pneumocystis murina (strain B123) (Mouse pneumocystis pneumonia agent) Msg Antibody Biotin Conjugated

Catalog Number: DZ41871-Biotin

### Description

Size	200 ul/vial
Form	Liquid
Purification	Immunogen affinity purified.
Contents	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% NaN <sub>3</sub> .
Conjugate	Biotin
Suggested Dilutions	Western blot, Optimal dilutions should be determined by end users. Immunohistochemistry (Paraffin-embedded Section), Optimal dilutions should be determined by end users. ELISA, Optimal dilutions should be determined by end users.

### Storage

At -20°C for one year from date of receipt. Avoid repeated freezing and thawing.

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Pneumocystis murina (strain B123) (Mouse pneumocystis pneumonia agent) Msg Antibody - Biotin

For Research Use Only. Not for use in diagnostic procedures.