

Anti-Pneumocystis murina (strain B123) (Mouse pneumocystis pneumonia agent) Msg Antibody Fluoro594 Conjugated

Catalog Number: DZ41871-Fluoro594

Description

Size	200 ul/vial
Form	Liquid
Purification	Immunogen affinity purified.
Contents	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% NaN ₃ .
Conjugate	Fluoro594 Excitation Wavelength: 593 nm Emission Wavelength: 618 nm
Suggested Dilutions	Flow Cytometry, Optimal dilutions should be determined by end users.

Storage

At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Pneumocystis murina (strain B123) (Mouse pneumocystis pneumonia agent) Msg Antibody - Fluoro594

For Research Use Only. Not for use in diagnostic procedures.