

Anti-Pneumocystis murina (strain B123) (Mouse pneumocystis pneumonia agent) Msg Antibody

Catalog Number: DZ41871-carrier-free

Overview

Product Name	Anti-Pneumocystis murina (strain B123) (Mouse pneumocystis pneumonia agent) Msg Antibody
Reactive Species	Pneumocystis murina
Description	Boster Bio Anti-Pneumocystis murina (strain B123) (Mouse pneumocystis pneumonia agent) Msg Antibody catalog # DZ41871. This antibody reacts with Pneumocystis murina.
Application	IF, IHC, WB
Clonality	Polyclonal
Formulation	Each vial contains 20mM PBS, 50% glycerol, 0.02% NaN ₃ .
Storage Instructions	At -20°C for one year, at 4°C for one month. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	M7NNF9

Technical Details

Form	Liquid
Concentration	500 ug/ml
Purification	Immunogen affinity purified.
Suggested Dilutions	Western blot, 0.25-0.5 ug/ml Immunohistochemistry (Paraffin-embedded Section), 2-5 ug/ml Immunofluorescence, 5 ug/ml

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Pneumocystis murina (strain B123) (Mouse pneumocystis pneumonia agent) Msg Antibody

For Research Use Only. Not for use in diagnostic procedures.