

## Anti-Danio rerio (Zebrafish) trh Antibody

Catalog Number: DZ41880

### Overview

Product Name	Anti-Danio rerio (Zebrafish) trh Antibody
Reactive Species	Zebrafish
Description	Boster Bio Anti-Danio rerio (Zebrafish) trh Antibody catalog # DZ41880. This antibody reacts with Zebrafish.
Application	IF
Clonality	Polyclonal
Formulation	Each vial contains 20mM PBS, 50% glycerol, 0.02% NaN3.
Storage Instructions	At -20°C for one year, at 4°C for one month. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	Q5EDF9

### Technical Details

Form	Liquid
Concentration	500 ug/ml
Purification	Immunogen affinity purified.
Suggested Dilutions	Immunofluorescence, 5 ug/ml

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Danio rerio (Zebrafish) trh Antibody

For Research Use Only. Not for use in diagnostic procedures.