

Anti-Plasminogen PLG Antibody Peroxidase Conjugated

Catalog Number: H00702

About PLG

Plasmin is released as a zymogen called plasminogen (PLG) from the liver into the systemic circulation. Two major glycoforms of plasminogen are present in humans. Type II plasminogen is preferentially recruited to the cell surface over the type I glycoform. Conversely, type I plasminogen appears more readily recruited to blood clots. In circulation, plasminogen adopts a closed, activation resistant conformation. Upon binding to clots, or to the cell surface, plasminogen adopts an open form that can be converted into active plasmin by a variety of enzymes, including tissue plasminogen activator (tPA), urokinase plasminogen activator (uPA), kallikrein, and factor XII (Hageman factor). Fibrin is a cofactor for plasminogen activation by tissue plasminogen activator. Urokinase plasminogen activator receptor (uPAR) is a cofactor for plasminogen activation by urokinase plasminogen activator. The conversion of plasminogen to plasmin involves the cleavage of the peptide bond between Arg-561 and Val-562. Deficiency in plasmin may lead to thrombosis, as clots are not degraded adequately. Plasminogen deficiency in mice leads to defective liver repair, defective wound healing, reproductive abnormalities. In humans, a rare disorder called plasminogen deficiency type I is caused by mutations of the PLG gene and is often manifested by ligneous conjunctivitis.


Overview

Product Name	Anti-Plasminogen PLG Antibody Peroxidase Conjugated
Reactive Species	Human
Description	Boster Bio Anti-Plasminogen PLG Antibody Peroxidase Conjugated (Catalog # H00702). Tested in ELISA, WB applications. This antibody reacts with Human.
Conjugate	HRP
Application	ELISA, IHC, WB
Clonality	Polyclonal
Formulation	0.02 M Potassium Phosphate, 0.15 M Sodium Chloride, pH 7.2, antibody formulated with stabilizing protein (10 mg/mL) - Immunoglobulin and Protease free, 0.01% (w/v) Gentamicin Sulfate. Do NOT add Sodium Azide! This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store vial at 4°C prior to restoration. For extended storage aliquot contents and freeze at -20°C or below. Avoid cycles of freezing and thawing. Centrifuge product if not completely clear after standing at room temperature. This product is stable for several weeks at 4°C as an undiluted liquid. Dilute only prior to immediate use. Expiration date is one (1) year from date of opening. (Ship at ambient temperature.)
Host	Goat
Uniprot ID	P00747

Technical Details

Immunogen	Plasminogen [Human Plasma]
Predicted Reactive Species	Bovine, Canine, Equine, Guinea Pig, Pig, Rabbit, Zebrafish
Isotype	IgG
Form	Lyophilized
Concentration	1.0 mg/mL by UV absorbance at 280 nm
Purification	Anti-Plasminogen is an IgG fraction antibody purified from monospecific antiserum by a multi-step process which includes delipidation, salt fractionation and ion exchange chromatography followed by extensive dialysis against the buffer stated above. Assay by immunoelectrophoresis resulted in a single precipitin arc against anti-Peroxidase, anti-Goat Serum as well as purified and partially purified Plasminogen [Human Plasma]. Cross-reactivity against Plasminogen from other sources is unknown.
Suggested Dilutions	<p>ELISA: 1:2,000 - 1:10,000 IHC: 1:250 - 1:1,000 WB: 1:500 - 1:2,500</p> <p>Anti-Plasminogen has been tested by ELISA and western blot. This product is assayed against 1.0 µg of Plasminogen [Human Plasma] in a standard capture ELISA using ABTS (2,2'-azino-bis-[3-ethylbenthiazoline-6-sulfonic acid]) code # ABTS-100 as a substrate for 30 minutes at room temperature. A working dilution of 1:1,000 to 1:4,000 of the reconstitution concentration is suggested for this product.</p>

Anti-Plasminogen PLG Antibody Peroxidase Conjugated (H00702) Images



Boster Goat anti Plasminogen antibody (200-101-208 lot 6571) was used to detect Plasminogen under reducing (R) and non-reducing (NR) conditions. Reduced samples of purified target proteins contained 4% BME and were boiled for 5 minutes. Samples of ~1ug of protein per lane were run by SDS-PAGE. Protein was transferred to nitrocellulose and probed with 1:3000 dilution of primary antibody. Detection shown was using Dylight 649 conjugated Donkey anti goat 1 hr RT. Images were collected using the BioRad VersaDoc System.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Plasminogen PLG Antibody Peroxidase Conjugated
For Research Use Only. Not for use in diagnostic procedures.