

Anti-beta Actin Mouse Monoclonal Antibody(HRP Conjugated)

Catalog Number: H01263


Overview

Product Name	Anti-beta Actin Mouse Monoclonal Antibody(HRP Conjugated)
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-beta Actin Mouse Monoclonal Antibody(HRP Conjugated) catalog # H01263. Tested in WB application. This antibody reacts with Human, Mouse, Rat.
Application	WB
Clonality	Monoclonal 17A12
Formulation	Rabbit IgG in stabilizing components, phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. *This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Mouse
Uniprot ID	P60709


Technical Details

Immunogen	A synthesized peptide derived from human beta Actin
Isotype	IgG2b
Form	Liquid
Concentration	0.5mg/ml
Purification	Affinity-chromatography
Suggested Dilutions	WB 1:500-2000


Anti-beta Actin Mouse Monoclonal Antibody(HRP Conjugated) (H01263) Images




Western blot analysis of beta Actin expression in (1) HeLa cell lysate, (2) Human liver lysate, (3) 3T3 cell lysate, (4) Mouse brain lysate, (5) PC12 cell lysate, (6) Rat heart lysate.



All lanes use the Antibody at 1:1W dilution for 1 hour at room temperature.



All lanes use the Antibody at 1:1W dilution for 1 hour at room temperature.



All lanes use the Antibody at 1:1W dilution for 1 hour at room temperature.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-beta Actin Mouse Monoclonal Antibody(HRP Conjugated)

For Research Use Only. Not for use in diagnostic procedures.