

AARSD1 (PTGES3L-AARSD1) (NM_025267) Human Over-expression Lysate

Catalog Number: LS100885850

Transient overexpression lysate of alanyl-tRNA synthetase domain containing 1 (AARSD1), transcript variant 2

Overview

Product Name	AARSD1 (PTGES3L-AARSD1) (NM_025267) Human Over-expression Lysate
Reactive Species	Human
Description	Transient overexpression lysate of alanyl-tRNA synthetase domain containing 1 (AARSD1), transcript variant 2
Expression Host	HEK293T
Tag	C-Myc/DDK
Detection Antibodies	Clone OTI4C5, Anti-DDK (FLAG) monoclonal antibody (M30971)
Components	For 100 µg, it contains two vials: 1 vial of 100 µg gene specific transient over-expression cell lysate in RIPA buffer 1 vial of 100 µg whole HEK293T cell lysate in RIPA buffer For 20 µg, it contains one vial: 1 vial of 20 ug lyophilized gene specific transient over-expression cell lysate
Storage Instructions	The lysate is shipped with dry ice. Upon receiving, store the sample at -80°C. Also after dilution, the protein sample should be aliquoted and stored at -80°C for long term storage. Avoid repeated freeze-thaw cycles. Lysate samples can be diluted with 2xSDS Sample Buffer provided. Lysate samples are stable for 12 months from the date of receipt when stored at -80°C.
Preparation	HEK293T cells in 10-cm dishes were transiently transfected with Transfection Reagent and 5ug TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs before collection. The cells were lysed in modified RIPA buffer (25mM Tris-HCl pH7.6, 150mM

NaCl, 1% NP-40, 1mM EDTA, 1xProteinase inhibitor cocktail mix, 1mM PMSF and 1mM Na₃VO₄), and then centrifuged to clarify the lysate. Protein concentration was measured by BCA kit. Cell lysates were aliquoted and stored at -20°C before shipping.

[AARSD1 \(PTGES3L-AARSD1\) \(NM_025267\) Human Over-expression Lysate](#)

For Research Use Only. Not for use in diagnostic procedures.