



Tissue Protein Lysate, Stomach

Catalog number: LS781000

This package insert must be read in its entirety before using this product.

For research use only. Not for use in diagnostic procedures.

Tissue Protein Lysate, Stomach

Catalog Number: LS781000, **Storage:** The lysate is shipped with dry ice. Store the tissue lysate at -80°C. These products are stable for one year from date of shipping when stored at -80°C without repeated freeze/thaws.

Overview

Product Name	Tissue Protein Lysate, Stomach
Sku	LS781000
Size	150 ug
Disease State	Normal
Diagnosis Category	Normal
Tissue	Stomach
Frozen Sample ID	FR00034396
Case ID	CI0000013144
Prep	3P
Vial ID	PL00002716
Tissue of Origin	Stomach
Site of Finding	Stomach
Appearance	Normal
Sample Diagnosis (from Pathology Verification)	Within normal limits
Normal	100%
Lesional	0%
Tumor	0%
Tumor Hypo/Acellular Stroma	0%
Tumor Hypercellular Stroma	0%
Necrosis	0%
Pathology Verification notes from H&E review	Non Tumor Structures: 90% Mucosa, 5% Submucosa, 5% Muscle
CASE Diagnosis (from medical center path report)	Adenocarcinoma of pancreas, ductal
Age	80
Gender	Female
Race	White or Caucasian
Tumor Grade	Not applicable
Minimum Stage Grouping	Not applicable
T (Extent of Primary Tumor)	Not applicable
N (Lymph node metastasis)	Not applicable
M (Distant metastasis)	Not applicable
Storage & Handling	The lysate is shipped with dry ice. Store the tissue lysate at -80°C. These products are stable for one year from date of shipping when stored at -80°C without repeated freeze/thaws.

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



For Research Use Only. Not for use in diagnostic procedures.