



Tissue Protein Lysate, Lung

Catalog number: LS781020

This package insert must be read in its entirety before using this product.

For research use only. Not for use in diagnostic procedures.

Tissue Protein Lysate, Lung

Catalog Number: LS781020, **Storage:** The lysate is shipped with dry ice. Store the tissue lysate at -80°C. These products are stable for one year from date of shipping when stored at -80°C without repeated freeze/thaws.

Overview

Product Name	Tissue Protein Lysate, Lung
SKU	LS781020
Size	150 ug
Disease State	Other Disease
Diagnosis Category	Fibrosis
Tissue	Lung
Frozen Sample ID	FR00034F65
Case ID	CI0000014888
Prep	2P
Vial ID	PL00002275
Tissue of Origin	Lung
Site of Finding	Lung
Appearance	Lesional
Sample Diagnosis (from Pathology Verification)	Fibrosis
Normal	0%
Lesional	100%
Tumor	0%
Tumor Hypo/Acellular Stroma	0%
Tumor Hypercellular Stroma	0%
Necrosis	0%
Pathology Verification notes from H&E review	Non Tumor Structures: 100% No normal architecture
CASE Diagnosis (from medical center path report)	Carcinoma of lung, squamous cell
Age	72
Gender	Male
Race	White or Caucasian
Tumor Grade	Not applicable
Minimum Stage Grouping	Not applicable
T (Extent of Primary Tumor)	Not applicable
N (Lymph node metastasis)	Not applicable
M (Distant metastasis)	Not applicable
Storage & Handling	The lysate is shipped with dry ice. Store the tissue lysate at -80°C. These products are stable for one year from date of shipping when stored at -80°C without repeated freeze/thaws.

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Tissue Protein Lysate, Lung

For Research Use Only. Not for use in diagnostic procedures.