



Tissue Protein Lysate, Colon: right

Catalog number: LS781049

This package insert must be read in its entirety before using this product.

For research use only. Not for use in diagnostic procedures.

Tissue Protein Lysate, Colon: right

Catalog Number: LS781049, **Storage:** The lysate is shipped with dry ice. Store the tissue lysate at -80°C. These products are stable for one year from date of shipping when stored at -80°C without repeated freeze/thaws.

Overview

Product Name	Tissue Protein Lysate, Colon: right
Sku	LS781049
Size	150 ug
Disease State	Normal
Diagnosis Category	Normal
Tissue	Colon
Frozen Sample ID	FR5B343162
Case ID	CI0000018187
Prep	2P
Vial ID	PL000020DE
Tissue of Origin	Colon: right
Site of Finding	Colon: right
Appearance	Normal
Sample Diagnosis (from Pathology Verification)	Within normal limits
Normal	100%
Lesional	0%
Tumor	0%
Tumor Hypo/Acellular Stroma	0%
Tumor Hypercellular Stroma	0%
Necrosis	0%
Pathology Verification notes from H&E review	Non Tumor Structures: 80% Mucosa, 10% Submucosa, 10% Muscle; Other Features/Comments: Lymphoid nodules in submucosa
CASE Diagnosis (from medical center path report)	Crohns disease
Age	70
Gender	Female
Race	White or Caucasian
Tumor Grade	Not applicable
Minimum Stage Grouping	Not applicable
T (Extent of Primary Tumor)	Not applicable
N (Lymph node metastasis)	Not applicable
M (Distant metastasis)	Not applicable
Storage & Handling	The lysate is shipped with dry ice. Store the tissue lysate at -80°C. These products are stable for one year from date of shipping when stored at -80°C without repeated freeze/thaws.

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



[Tissue Protein Lysate, Colon: right](#)

For Research Use Only. Not for use in diagnostic procedures.