



Tissue Protein Lysate, Lung: right lower lobe

Catalog number: LS781053

This package insert must be read in its entirety before using this product.

For research use only. Not for use in diagnostic procedures.

Tissue Protein Lysate, Lung: right lower lobe

Catalog Number: LS781053, **Storage:** The lysate is shipped with dry ice. Store the tissue lysate at -80°C. These products are stable for one year from date of shipping when stored at -80°C without repeated freeze/thaws.

Overview

Product Name	Tissue Protein Lysate, Lung: right lower lobe
Sku	LS781053
Size	150 ug
Disease State	Cancer
Diagnosis Category	Lung Cancer
Tissue	Lung
Frozen Sample ID	FR5B33DC9F
Case ID	CI0000018513
Prep	1P
Vial ID	PL00002671
Tissue of Origin	Lung: right lower lobe
Site of Finding	Lung: right lower lobe
Appearance	Tumor
Sample Diagnosis (from Pathology Verification)	Adenocarcinoma of lung, papillary
Normal	0%
Lesional	0%
Tumor	70%
Tumor Hypo/Acellular Stroma	0%
Tumor Hypercellular Stroma	30%
Necrosis	0%
CASE Diagnosis (from medical center path report)	Adenocarcinoma of lung, papillary
Age	66
Gender	Male
Race	White or Caucasian
Tumor Grade	AJCC G2: Moderately differentiated
Minimum Stage Grouping	IIIA
T (Extent of Primary Tumor)	T2
N (Lymph node metastasis)	N2
M (Distant metastasis)	MX
Storage & Handling	The lysate is shipped with dry ice. Store the tissue lysate at -80°C. These products are stable for one year from date of shipping when stored at -80°C without repeated freeze/thaws.

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon giftcard! Your



reviews help your fellow scientists make the right decisions. Thank you for your contribution.

[Tissue Protein Lysate, Lung: right lower lobe](#)

For Research Use Only. Not for use in diagnostic procedures.