

## Anti-ERBB2 / HER2 / CD340 Reference Antibody (disitamab)

Catalog Number: M00010-11


### Overview

Product Name	Anti-ERBB2 / HER2 / CD340 Reference Antibody (disitamab)
Reactive Species	Cynomolgus monkey, Human, Mouse, Rat
Description	Boster Bio Anti-ERBB2 / HER2 / CD340 Reference Antibody (disitamab) (Catalog # M00010-11). Tested in Flow Cytometry, ELISA, FTA. This antibody reacts with Rat, Cynomolgus monkey, Human, Mouse. Endotoxin: < 0.033EU/ug,determined by LAL method. Expression system: CHO Cell
Application	ELISA, Flow Cytometry, Functional Assay, Kinetics
Clonality	Monoclonal
Formulation	Purified monoclonal antibody supplied in PBS, pH6.0, without preservative.This antibody is purified through a protein A column.
Storage Instructions	Maintain refrigerated at 2-8°C for up to 2 weeks. For long-term storage, aliquot and store at -20°C to avoid repeated freeze-thaw cycles.
Uniprot ID	P04626


### Technical Details

Isotype	IgG1
Form	Liquid


## Anti-ERBB2 / HER2 / CD340 Reference Antibody (disitamab) (M00010-11) Images



Anti-ERBB2/HER2/CD340 Reference Antibody (disitamab) on SDS-PAGE under reducing (R) condition. The gel was stained with Coomassie Blue. The purity of the protein is greater than 95%



The purity of Anti-ERBB2/HER2/CD340 Reference Antibody (disitamab) is more than 98.73%



Immobilized human ERBB2/HER2/CD340

### Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-ERBB2 / HER2 / CD340 Reference Antibody (disitamab)

For Research Use Only. Not for use in diagnostic procedures.