

Anti-TERT Rabbit Monoclonal Antibody

Catalog Number: M00012-2


Overview

Product Name	Anti-TERT Rabbit Monoclonal Antibody
Reactive Species	Human
Description	Boster Bio Anti-TERT Rabbit Monoclonal Antibody catalog # M00012-2. Tested in WB, IP applications. This antibody reacts with Human.
Application	IP, WB
Clonality	Monoclonal 17T17
Formulation	Rabbit IgG in stabilizing components, phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. *This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	O14746

Technical Details

Immunogen	A synthesized peptide derived from human TERT
Isotype	IgG
Form	Liquid
Concentration	0.5mg/ml
Purification	Affinity-chromatography
Suggested Dilutions	WB 1:500-2000 IP 1:50-100

Anti-TERT Rabbit Monoclonal Antibody (M00012-2) Images



Western blot analysis of TERT expression in HeLa cells lysate.

4 Publications Citing This Product

1. PubMed ID: 10.3892/or.19.2.575, Long-term effects of short hairpin RNA-targeted human telomerase reverse transcriptase on suppression of SGC-7901 cell proliferation by inhibition of telomerase activity
2. PubMed ID: 10.4161/cbt.8.10.8117, Adenoviral mediated transduction of adenoid cystic carcinoma by human TRAIL gene driven with hTERT tumor specific promoter induces apoptosis.
3. PubMed ID: 10.1007/s12250-008-2914-2, Expression of human papillomavirus in mammary carcinoma and its possible mechanism in carcinogenesis

Visit bosterbio.com/anti-tert-rabbit-monoclonal-antibody-m00012-2-boster.html to see all 4 publications.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-TERT Rabbit Monoclonal Antibody

For Research Use Only. Not for use in diagnostic procedures.