

## Anti-Akt 1/2/3 Rabbit Monoclonal Antibody

Catalog Number: M00024-10


### Overview

|                      |  |
|----------------------|--|
| Product Name         | Anti-Akt 1/2/3 Rabbit Monoclonal Antibody  |
| Reactive Species     | Human, Mouse   |
| Description          | This Recombinant Rabbit Monoclonal Antibody detects endogenous levels of Akt 1/2/3 protein. Validated for research with WB, IP,ELISA applications. Anti Akt 1/2/3 is reactive for Human, Mouse, (predicted: Rat) samples   |
| Application          | ELISA, IP, WB  |
| Clonality            | Monoclonal   |
| Formulation          | PBS, 50% glycerol, 0.05% Proclin 300, 0.05% stabilizing protein<br>This antibody is supplied in a stabilized formulation.<br>Compatibility with conjugation reactions depends on the chemistry of the conjugation method used.<br>For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required. |
| Storage Instructions | Store at -20°C/1 year  |
| Host                 | Rabbit   |
| Uniprot ID           | P31749/P31751/Q9Y243   |

### Technical Details

|                            |   |
|----------------------------|---|
| Predicted Reactive Species | Rat                                     |
| Isotype                    | IgG, Kappa                              |
| Purification               | Protein A                               |
| Suggested Dilutions        | WB 1:500-1:2000<br>ELISA 1:5000-1:20000 |

## Anti-Akt 1/2/3 Rabbit Monoclonal Antibody (M00024-10) Images



Various whole cell lysates were separated by 10% SDS-PAGE, and the membrane was blotted with anti-AKT antibody. The HRP-conjugated Goat anti-Rabbit IgG(H + L) antibody was used to detect the antibody. Lane 1:mouse testis Lane 2:HeLa Predicted band size: 60kDa Observed band size: 60kDa

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Akt 1/2/3 Rabbit Monoclonal Antibody

For Research Use Only. Not for use in diagnostic procedures.