

Anti-AKT (PH domain) AKT1 Rabbit Monoclonal Antibody, Clone#RM316

Catalog Number: M00024-4


Overview

Product Name	Anti-AKT (PH domain) AKT1 Rabbit Monoclonal Antibody, Clone#RM316
Reactive Species	Human
Description	Boster Bio Anti-AKT (PH domain) AKT1 Rabbit Monoclonal Antibody, Clone#RM316 (Catalog # M00024-4). Tested in IHC, WB applications. This antibody reacts with Human.
Application	IHC, WB
Clonality	Monoclonal RM316
Formulation	50% Glycerol/PBS with 1% stabilizing protein and 0.09% sodium azide This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P31749


Technical Details

Immunogen	A peptide corresponding to the PH domain of human AKT
Cross Reactivity	This antibody reacts to human AKT (PH domain). It may also react to mouse or rat AKT, as predicted by immunogen homology.
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5-1mg/ml, actual concentration vary by lot. Use suggested dilution ratio to decide dilution procedure.
Purification	Protein A affinity purified from an animal origin-free culture supernatant
Suggested Dilutions	Immunohistochemistry (IHC): 1:200-1:500 dilution WB: 1:5000-1:10000 dilution.

Anti-AKT (PH domain) AKT1 Rabbit Monoclonal Antibody, Clone#RM316 (M00024-4) Images



Western Blotting result Western Blot of 293 cells lysate using anti-AKT (PH domain) rabbit monoclonal antibody (Clone RM316) at a 1:10000 dilution.



IHC result Immunohistochemical staining of formalin fixed and paraffin embedded human breast cancer tissue section using anti-AKT (PH domain) rabbit monoclonal antibody (Clone RM316) at a 1:500 dilution.

2 Publications Citing This Product

1. PubMed ID: 33382670, Chen H, Sheng H, Zhao Y, Zhu G. Piperine Inhibits Cell Proliferation and Induces Apoptosis of Human Gastric Cancer Cells by Downregulating Phosphatidylinositol 3-Kinase (PI3K)/Akt Pathway. Med Sci Monit. 2020 Dec 31;26:e928403. doi:10.12659/MSM.928403. PMID:33382

2. PubMed ID: 18694621, A small interfering RNA targeting osteopontin as gastric cancer therapeutics

Visit [bosterbio.com/anti-akt-ph-domain-rabbit-monoclonal-antibody-clone-rm316-m00024-4-boster.html](https://www.bosterbio.com/anti-akt-ph-domain-rabbit-monoclonal-antibody-clone-rm316-m00024-4-boster.html) to see all 2 publications.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-AKT (PH domain) AKT1 Rabbit Monoclonal Antibody, Clone#RM316

For Research Use Only. Not for use in diagnostic procedures.