

Anti-Myc-Tag Mouse Monoclonal Antibody

Catalog Number: M00026-5


Overview

Product Name	Anti-Myc-Tag Mouse Monoclonal Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-Myc-Tag Mouse Monoclonal Antibody catalog # M00026-5. Tested in WB, ICC/IF, IP applications. This antibody reacts with Human, Mouse, Rat.
Application	IP, IF, ICC, WB
Clonality	Monoclonal 17M02
Formulation	Mouse IgG in stabilizing components, phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. *This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Mouse
Uniprot ID	N/A

Technical Details

Immunogen	A synthesized peptide
Isotype	IgG1
Form	Liquid
Concentration	0.5mg/ml
Purification	Affinity-chromatography
Suggested Dilutions	WB 1:1000-1:5000 ICC/IF 1:100-1:200 IP 1:100

Anti-Myc-Tag Mouse Monoclonal Antibody (M00026-5) Images



Western blot analysis of Myc-tagged protein expression in (1) 293T cell lysate; (2) 293T cell transfected with Myc-tagged protein lysate.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Myc-Tag Mouse Monoclonal Antibody

For Research Use Only. Not for use in diagnostic procedures.