

Anti-VEGF Reference Antibody (BioMab patent anti-VEGF)

Catalog Number: M00045-3

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.




Overview

Product Name	Anti-VEGF Reference Antibody (BioMab patent anti-VEGF)
Reactive Species	Human, Mouse
Description	Boster Bio Anti-VEGF Reference Antibody (BioMab patent anti-VEGF) (Catalog # M00045-3). Tested in Flow Cytometry, ELISA, FTA. This antibody reacts with Human, Mouse. Endotoxin: < 0.086EU/ug,determined by LAL method. Expression system: CHO Cell
Application	ELISA, Flow Cytometry, Functional Assay, Kinetics
Clonality	Monoclonal
Formulation	Purified monoclonal antibody supplied in PBS, pH6.0, without preservative.This antibody is purified through a protein A column.
Storage Instructions	Maintain refrigerated at 2-8°C for up to 2 weeks. For long-term storage, aliquot and store at -20°C to avoid repeated freeze-thaw cycles.
Uniprot ID	P15692


Technical Details

Isotype	IgG1
Form	Liquid


Anti-VEGF Reference Antibody (BioMab patent anti-VEGF) (M00045-3) Images




Anti-VEGF Reference Antibody (BioMab patent anti-VEGF) on SDS-PAGE under reducing (R) condition. The gel was stained with Coomassie Blue. The purity of the protein is greater than 95%



The purity of Anti-VEGF Reference Antibody (BioMab patent anti-VEGF) is more than 97.8%



Immobilized human VEGF165 His at 2 µg/mL can bind Anti-VEGF Reference Antibody (BioMab patent anti-VEGF)



Bevacizumab inhibited the tumor growth of COLO205 on balb/c nude mice. The result showed significant anti-tumor effects

Anti-VEGF Reference Antibody (BioMab patent anti-VEGF)

For Research Use Only. Not for use in diagnostic procedures.