

Anti-Phospho-MDM2 (Tyr394) Antibody

Catalog Number: M00054Y394

Overview

Product Name	Anti-Phospho-MDM2 (Tyr394) Antibody
Reactive Species	Human, Mouse
Description	Boster Bio Anti-Phospho-MDM2 (Tyr394) Antibody catalog # M00054Y394. Tested in ELISA, WB applications. This antibody reacts with Human, Mouse.
Application	ELISA, WB
Clonality	Polyclonal
Formulation	Liquid in PBS containing 50% glycerol, 0.5% stabilizing protein and 0.02% sodium azide. This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q00987

Technical Details

Immunogen	Synthesized phospho-peptide around the phosphorylation site of human MDM2 (Phospho-Tyr394)
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Suggested Dilutions	WB 1:500-2000 ELISA 1:10000-20000

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your



reviews help your fellow scientists make the right decisions. Thank you for your contribution.

Anti-Phospho-MDM2 (Tyr394) Antibody

For Research Use Only. Not for use in diagnostic procedures.