

Anti-IGF1R / CD221 Reference Antibody (Immunomedics patent anti-IGF-1R)

Catalog Number: M00070-11


Overview

Product Name	Anti-IGF1R / CD221 Reference Antibody (Immunomedics patent anti-IGF-1R)
Reactive Species	Human
Description	Boster Bio Anti-IGF1R / CD221 Reference Antibody (Immunomedics patent anti-IGF-1R) (Catalog # M00070-11). Tested in Flow Cytometry, ELISA, FTA. This antibody reacts with Human. Endotoxin: < 0.121EU/ug,determined by LAL method. Expression system: CHO Cell
Application	ELISA, Flow Cytometry, Functional Assay, Kinetics
Clonality	Monoclonal
Formulation	Purified monoclonal antibody supplied in PBS, pH6.0, without preservative.This antibody is purified through a protein A column.
Storage Instructions	Maintain refrigerated at 2–8°C for up to 2 weeks. For long-term storage, aliquot and store at -20°C to avoid repeated freeze-thaw cycles.
Uniprot ID	P08069


Technical Details

Isotype	IgG1
Form	Liquid

Anti-IGF1R / CD221 Reference Antibody (Immunomedics patent anti-IGF-1R) (M00070-11) Images



Anti-IGF1R/CD221 Reference Antibody (Immunomedics patent anti-IGF-1R) on SDS-PAGE under reducing (R) condition. The gel was stained with Coomassie Blue. The purity of the protein is greater than 95%



The purity of Anti-IGF1R/CD221 Reference Antibody (Immunomedics patent anti-IGF-1R) is more than 95%

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-IGF1R / CD221 Reference Antibody (Immunomedics patent anti-IGF-1R)

For Research Use Only. Not for use in diagnostic procedures.