

Anti-B Raf Rabbit Monoclonal Antibody

Catalog Number: M00075


Overview

Product Name	Anti-B Raf Rabbit Monoclonal Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-B Raf Rabbit Monoclonal Antibody catalog # M00075. Tested in WB, IHC, IP, Flow Cytometry applications. This antibody reacts with Human, Mouse, Rat.
Application	Flow Cytometry, IP, IHC, WB
Clonality	Monoclonal CCG-2
Formulation	Rabbit IgG in stabilizing components, phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. *This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P15056


Technical Details

Immunogen	A synthesized peptide derived from human B Raf
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5mg/ml
Purification	Affinity-chromatography
Suggested Dilutions	WB 1:500-2000 IHC 1:50-200 IP 1:20 FC 1:20


Anti-B Raf Rabbit Monoclonal Antibody (M00075) Images



All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.



Western blot analysis of B Raf expression in (1)Human brain lysate (2)HepG2 cell lysate.



Immunohistochemical analysis of paraffin-embedded human ovarian cancer, using B Raf Antibody.

1 Publications Citing This Product

1. PubMed ID: -, Huanyue Cui,Xueying Liu,Jin Zhang,Ke Zhang,Dahong Yao,Shi Dong,Shushu Feng,Lu Yang, Yuyao Li,Hangyu Wang,Jian Huang,Jinhui Wang,"Rhodiola rosea L. Attenuates Cigarette Smoke and Lipopolysaccharide-Induced COPD in Rats via Inflammation Inhibition and Antioxidant and Antifibrosis Pathways",Evidence-Based Complementary and Alternative Medicine,vol. 2021,Article ID 6103158,18 pages,2021.<https://doi.org/10.1155/2021/6103158>

Visit [bosterbio.com/anti-b-raf-rabbit-monoclonal-antibody-m00075-boster.html](https://www.bosterbio.com/anti-b-raf-rabbit-monoclonal-antibody-m00075-boster.html) to see all 1 publications.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-B Raf Rabbit Monoclonal Antibody

For Research Use Only. Not for use in diagnostic procedures.