

Anti-LDL Receptor Rabbit Monoclonal Antibody

Catalog Number: M00076


Overview

Product Name	Anti-LDL Receptor Rabbit Monoclonal Antibody
Reactive Species	Human, Mouse
Description	Boster Bio Anti-LDL Receptor Rabbit Monoclonal Antibody catalog # M00076. Tested in WB, IHC applications. This antibody reacts with Human, Mouse.
Application	IHC, WB
Clonality	Monoclonal BFC-12
Formulation	Rabbit IgG in stabilizing components, phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. *This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P01130


Technical Details

Immunogen	A synthesized peptide derived from human LDL Receptor
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5mg/ml
Purification	Affinity-chromatography
Suggested Dilutions	WB 1:500-2000 IHC 1:50-200

Anti-LDL Receptor Rabbit Monoclonal Antibody (M00076) Images



Western blot analysis of LDLR expression in HepG2 cell lysate.



Immunohistochemical analysis of paraffin-embedded human liver carcinoma, using LDL Receptor Antibody .

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-LDL Receptor Rabbit Monoclonal Antibody

For Research Use Only. Not for use in diagnostic procedures.