

Anti-COX-2 PTGS2 Rabbit Monoclonal Antibody, Clone#RM348

Catalog Number: M00084-2


Overview

Product Name	Anti-COX-2 PTGS2 Rabbit Monoclonal Antibody, Clone#RM348
Reactive Species	Human
Description	Boster Bio Anti-COX-2 PTGS2 Rabbit Monoclonal Antibody, Clone#RM348 (Catalog # M00084-2). Tested in IHC, WB applications. This antibody reacts with Human.
Application	IHC, WB
Clonality	Monoclonal RM348
Formulation	50% Glycerol/PBS with 1% stabilizing protein and 0.09% sodium azide This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P35354


Technical Details

Immunogen	A peptide corresponding to the C-terminus of human COX-2.
Cross Reactivity	This antibody reacts to human COX-2 (Prostaglandin G/H synthase 2).
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5-1mg/ml, actual concentration vary by lot. Use suggested dilution ratio to decide dilution procedure.
Purification	Protein A affinity purified from an animal origin-free culture supernatant
Suggested Dilutions	Immunohistochemistry (IHC): 1:500-1:1000 dilution WB: 1:1000-1:20000 dilution.

Anti-COX-2 PTGS2 Rabbit Monoclonal Antibody, Clone#RM348 (M00084-2) Images



Western Blotting result Western Blot of HeLa cell lysate using anti-COX-2 rabbit monoclonal antibody (Clone RM348) at a 1:20000 dilution.



IHC result Immunohistochemical staining of formalin fixed and paraffin embedded human lung cancer tissue section using anti-COX-2 rabbit monoclonal antibody (Clone RM348) at a 1:1000 dilution.

1 Publications Citing This Product

1. PubMed ID: 26617722, MicroRNA-101 regulates the viability and invasion of cervical cancer cells

Visit bosterbio.com/anti-cox-2-rabbit-monoclonal-antibody-clone-rm348-m00084-2-boster.html to see all 1 publications.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-COX-2 PTGS2 Rabbit Monoclonal Antibody, Clone#RM348

For Research Use Only. Not for use in diagnostic procedures.