

Anti-Rad51 Rabbit Monoclonal Antibody

Catalog Number: M00088


Overview

| | |
|----------------------|--|
| Product Name | Anti-Rad51 Rabbit Monoclonal Antibody |
| Reactive Species | Human, Mouse, Rat |
| Description | Boster Bio Anti-Rad51 Rabbit Monoclonal Antibody catalog # M00088. Tested in WB, Flow Cytometry, IP applications. This antibody reacts with Human, Mouse, Rat. |
| Application | Flow Cytometry, IP, WB |
| Clonality | Monoclonal CAE-18 |
| Formulation | Rabbit IgG in stabilizing components, phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. *This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required. |
| Storage Instructions | Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles. |
| Host | Rabbit |
| Uniprot ID | Q06609 |

Technical Details

| | |
|---------------------|--|
| Immunogen | A synthesized peptide derived from human Rad51 |
| Isotype | Rabbit IgG |
| Form | Liquid |
| Concentration | 0.5mg/ml |
| Purification | Affinity-chromatography |
| Suggested Dilutions | WB 1:500-1:2000 IP 1:50 FC 1:50 |

Anti-Rad51 Rabbit Monoclonal Antibody (M00088) Images



Western blot analysis of Rad51 expression in HeLa cell lysate.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Rad51 Rabbit Monoclonal Antibody

For Research Use Only. Not for use in diagnostic procedures.