

Anti-Transferrin Purified Azide Free TF Monoclonal Antibody

Catalog Number: M00094-1

About TF

Transferrin is a monomeric glycoprotein of approximately 77 kDa, which serves as an iron-transporter. In normal plasma, transferrin has a concentration of 25-50 μmol / liter, and is usually about one-third saturated with iron, thus providing a large buffering capacity in case of an acute increase in plasma iron levels. Cells take up transferrin-iron complexes (holotransferrin) using transferrin receptor dimers. Upon binding of holotransferrin, the receptor is internalized by clathrin-mediated endocytosis. Acidification of endosomes by vesicular membrane proton pumps leads to dissociation of iron ions, whereas transferrin (apotransferrin) remains associated with its receptor (CD71) and recycles to the cell surface, where apotransferrin is released upon exposure to normal pH. Internalization of labeled transferrin thus represents an useful approach to study endocytosis. Serum concentration rises in iron deficiency and pregnancy and falls in iron overload, infection and inflammatory conditions. Iron/transferrin complex is essential in haemoglobin synthesis and for certain types of cell division.

Overview


Product Name	Anti-Transferrin Purified Azide Free TF Monoclonal Antibody
Reactive Species	Human, Pig, Rabbit
Description	Boster Bio Anti-Transferrin Purified Azide Free TF Monoclonal Antibody (Catalog# M00094-1). Tested in WB, IHC-P, ICC, ELISA, RIA, Functional Assay application(s). This antibody reacts with Human, Pig, Rabbit.
Application	ELISA, IHC-P, ICC, WB, Functional Assay, RIA
Clonality	Monoclonal HTF-14
Formulation	Phosphate buffered saline (PBS), pH 7.4, azide-free
Storage Instructions	Store at 2-8°C. Do not freeze.
Host	Mouse
Uniprot ID	P02787

Technical Details


Immunogen	Purified porcine transferrin. The antibody HTF-14 recognizes an epitope located in the N-terminal domain of human serum transferrin, a 77 kDa single polypeptide chain glycoprotein (member of the iron binding family of proteins). It is synthesised in the liver and consists of two domains each having a high affinity reversible binding site for Fe ³⁺ .
Predicted Reactive Species	Primate
Isotype	Mouse IgG1
Form	Liquid

Concentration	1 mg/ml
Purification	Purified by protein-A affinity chromatography.
Suggested Dilutions	Functional application: The antibody HTF-14 blocks binding of transferrin to its receptor. Immunohistochemistry (paraffin sections): 10 ug/ml; positive tissue: placenta. Western blotting: non-reducing conditions, recommended dilution: 1-2 ug/ml.


Anti-Transferrin Purified Azide Free TF Monoclonal Antibody (M00094-1) Images



Human transferrin detected by the mouse monoclonal antibody HTF-14. 1. hTransferrin; 5 µg/well (red. con.) 2. hTransferrin; 3 µg/well (red. con.) 3. hTransferrin; 1 µg/well (red. con.) M Low Range marker (Bio-Rad) 4. hTransferrin; 1 µg/well (non-red. con.) 5. hTransferrin; 3 µg/well (non-red. con.) 6. hTransferrin; 5 µg/well (non-red. con.)



Human transferrin detected by the mouse monoclonal antibody HTF-14. 1. hTransferrin; 5 µg/well (red. con.) 2. hTransferrin; 3 µg/well (red. con.) 3. hTransferrin; 1 µg/well (red. con.) M Low Range marker (Bio-Rad) 4. hTransferrin; 1 µg/well (non-red. con.) 5. hTransferrin; 3 µg/well (non-red. con.) 6. hTransferrin; 5 µg/well (non-red. con.)



Immunohistochemistry staining of human placenta (paraffin sections) using anti-transferrin (HTF-14). Commercially tested by LifeSpan BioSciences.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Transferrin Purified Azide Free TF Monoclonal Antibody

For Research Use Only. Not for use in diagnostic procedures.