

## Anti-PDGFR-beta Rabbit Monoclonal Antibody

Catalog Number: M00096-6


### Overview

Product Name	Anti-PDGFR-beta Rabbit Monoclonal Antibody
Reactive Species	Human, Mouse, Rat
Description	This Recombinant Rabbit Monoclonal Antibody detects endogenous levels of PDGFR-beta protein. Validated for research with WB,IHC,IF,IP,ELISA applications. Anti PDGFR-beta is reactive for Human, Mouse, Rat samples
Application	ELISA, IP, IF, IHC, WB
Clonality	Monoclonal
Formulation	PBS, 50% glycerol, 0.05% Proclin 300, 0.05% stabilizing protein This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C/1 year
Host	Rabbit
Uniprot ID	P09619


### Technical Details

Isotype	IgG, Kappa
Purification	Protein A
Suggested Dilutions	IHC 1:200-1:1000 WB 1:2000-1:10000 IF 1:200-1:1000 ELISA 1:5000-1:20000 IP 1:50-1:200


## Anti-PDGFR-beta Rabbit Monoclonal Antibody (M00096-6) Images




Mouse spleen was stained with anti-PDGFR-beta rabbit antibody



Rat spleen was stained with anti-PDGFR-beta rabbit antibody



Human kidney was stained with anti-PDGFR-beta rabbit antibody



Various whole cell lysates were separated by 4-8% SDS-PAGE, and the membrane was blotted with anti-PDGFR-beta antibody. The HRP-conjugated Goat anti-Rabbit IgG(H + L) antibody was used to detect the antibody. Lane 1: C6 Lane 2: NIH-3T3 Lane 3: MG-63 Predicted band size: 124kDa Observed band size: 150-190kDa

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-PDGFR- $\beta$  Rabbit Monoclonal Antibody

For Research Use Only. Not for use in diagnostic procedures.