

Anti-FGFR1 Rabbit Monoclonal Antibody

Catalog Number: M00098


Overview

Product Name	Anti-FGFR1 Rabbit Monoclonal Antibody
Reactive Species	Human
Description	Boster Bio Anti-FGFR1 Rabbit Monoclonal Antibody catalog # M00098. Tested in WB, ICC/IF applications. This antibody reacts with Human.
Application	IF, ICC, WB
Clonality	Monoclonal EHC-6
Formulation	Rabbit IgG in stabilizing components, phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. *This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P11362


Technical Details

Immunogen	A synthesized peptide derived from human FGFR1
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5mg/ml
Purification	Affinity-chromatography
Suggested Dilutions	WB 1:500-2000 ICC/IF 1:50-200


Anti-FGFR1 Rabbit Monoclonal Antibody (M00098) Images




Immunohistochemical analysis of paraffin-embedded Human pituitary tumor, using the Antibody at 1:50 dilution.



Immunofluorescent analysis using the Antibody at 1:150 dilution.





Immunohistochemical analysis of paraffin-embedded Human thymoma, using the Antibody at 1:50 dilution.



Western blot analysis of FGFR1 expression in SH-SY5Y cell lysate.

Immunohistochemical analysis of paraffin-embedded Rat pancreas, using the Antibody at 1:50 dilution.



Immunofluorescent analysis of HeLa cells, using FGFR1 Antibody.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-FGFR1 Rabbit Monoclonal Antibody

For Research Use Only. Not for use in diagnostic procedures.