

Anti-IL-1b Reference Antibody (canakinumab)


Catalog Number: M00101-5

Overview


Product Name	Anti-IL-1b Reference Antibody (canakinumab)
Reactive Species	Human, Mouse
Description	Boster Bio Anti-IL-1b Reference Antibody (canakinumab) (Catalog # M00101-5). Tested in Flow Cytometry, ELISA, FTA. This antibody reacts with Human, Mouse. Endotoxin: < 0.007EU/ug,determined by LAL method. Expression system: CHO Cell
Application	ELISA, Flow Cytometry, Functional Assay, Kinetics
Clonality	Monoclonal
Storage Instructions	Maintain refrigerated at 2-8°C for up to 2 weeks. For long-term storage, aliquot and store at -20°C to avoid repeated freeze-thaw cycles.
Uniprot ID	P01584

Technical Details

Anti-IL-1b Reference Antibody (canakinumab) (M00101-5) Images



Anti-IL-1b Reference Antibody (canakinumab) on SDS-PAGE under reducing (R) condition. The gel was stained with Coomassie Blue. The purity of the protein is greater than 95%



The purity of Anti-IL-1b Reference Antibody (canakinumab) is more than 100%

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-IL-1b Reference Antibody (canakinumab)

For Research Use Only. Not for use in diagnostic procedures.