

Anti-IL-6Ra / CD126 Reference Antibody (vobarilizumab)

Catalog Number: M00102-7


Overview

Product Name	Anti-IL-6Ra / CD126 Reference Antibody (vobarilizumab)
Reactive Species	Cynomolgus monkey, Human
Description	Boster Bio Anti-IL-6Ra / CD126 Reference Antibody (vobarilizumab) (Catalog # M00102-7). Tested in Flow Cytometry, ELISA, FTA. This antibody reacts with Cynomolgus monkey, Human. Endotoxin: < 0.185EU/ug,determined by LAL method. Expression system: CHO Cell
Application	ELISA, Flow Cytometry, Functional Assay, Kinetics
Clonality	Monoclonal
Formulation	Purified monoclonal antibody supplied in PBS, pH6.0, without preservative.This antibody is purified through a protein A column.
Storage Instructions	Maintain refrigerated at 2-8°C for up to 2 weeks. For long-term storage, aliquot and store at -20°C to avoid repeated freeze-thaw cycles.
Uniprot ID	P08887


Technical Details

Isotype	IgG1
Form	Liquid


Anti-IL-6Ra / CD126 Reference Antibody (vobarilizumab) (M00102-7) Images




Anti-IL-6Ra/CD126 Reference Antibody (vobarilizumab) on SDS-PAGE under reducing (R) condition. The gel was stained with Coomassie Blue. The purity of the protein is greater than 96.9%



The purity of Anti-IL-6Ra/CD126 Reference Antibody (vobarilizumab) is more than 99.34%



Immobilized human IL 6R His at 2 µg/mL can bind Anti-IL-6Ra/CD126 Reference Antibody (vobarilizumab)



Human IL6R CHO cells were stained with Anti-IL-6Ra/CD126 Reference Antibody (vobarilizumab) and negative control protein respectively

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-IL-6Ra / CD126 Reference Antibody (vobarilizumab)

For Research Use Only. Not for use in diagnostic procedures.