

Anti-YAP Rabbit Monoclonal Antibody

Catalog Number: M00116-1

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.




Overview

Product Name	Anti-YAP Rabbit Monoclonal Antibody
Reactive Species	Human, Mouse, Rat
Description	This Recombinant Rabbit Monoclonal Antibody detects endogenous levels of YAP protein. Validated for research with WB,IHC,IF,IP,ELISA applications. Anti YAP is reactive for Human, Mouse, Rat samples
Application	ELISA, IP, IF, IHC, WB
Clonality	Monoclonal
Formulation	PBS, 50% glycerol, 0.05% Proclin 300, 0.05% stabilizing protein This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C/1 year
Host	Rabbit
Uniprot ID	P46937


Technical Details

Isotype	IgG, Kappa
Purification	Protein A
Suggested Dilutions	IHC 1:100-1:500 WB 1:500-1:2000 IF 1:200-1:1000 ELISA 1:5000-1:20000 IP 1:50-1:200


Anti-YAP Rabbit Monoclonal Antibody (M00116-1) Images




Mouse thyroid was stained with anti-YAP rabbit antibody




Rat thyroid was stained with anti-YAP rabbit antibody



Mouse skin was stained with anti-YAP rabbit antibody



Human breast carcinoma was stained with anti-YAP rabbit antibody



Various whole cell lysates were separated by 4-20% SDS-PAGE, and the membrane was blotted with anti-YAP antibody. The HRP-conjugated Goat anti-Rabbit IgG(H + L) antibody was used to detect the antibody. Lane 1: MCF7 Lane 2: 4T1 Lane 3: C6 Predicted band size: 55kDa Observed band size: 55-75kDa

Anti-YAP Rabbit Monoclonal Antibody

For Research Use Only. Not for use in diagnostic procedures.