

Anti-TLR2 / CD282 Reference Antibody (tomaralimab)

Catalog Number: M00131-2


Overview

Product Name	Anti-TLR2 / CD282 Reference Antibody (tomaralimab)
Reactive Species	Human, Mouse
Description	Boster Bio Anti-TLR2 / CD282 Reference Antibody (tomaralimab) (Catalog # M00131-2). Tested in Flow Cytometry, ELISA, FTA. This antibody reacts with Human, Mouse. Endotoxin: < 0.216EU/ug,determined by LAL method. Expression system: CHO Cell
Application	ELISA, Flow Cytometry, Functional Assay, Kinetics
Clonality	Monoclonal
Formulation	Purified monoclonal antibody supplied in PBS, pH6.0, without preservative.This antibody is purified through a protein A column.
Storage Instructions	Maintain refrigerated at 2-8°C for up to 2 weeks. For long-term storage, aliquot and store at -20°C to avoid repeated freeze-thaw cycles.
Uniprot ID	O60603


Technical Details

Isotype	IgG4
Form	Liquid

Anti-TLR2 / CD282 Reference Antibody (tomaralimab) (M00131-2) Images



Anti-TLR2/CD282 Reference Antibody (tomaralimab) on SDS-PAGE under reducing (R) condition. The gel was stained with Coomassie Blue. The purity of the protein is greater than 95%



The purity of Anti-TLR2/CD282 Reference Antibody (tomaralimab) is more than 100%

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-TLR2 / CD282 Reference Antibody (tomaralimab)

For Research Use Only. Not for use in diagnostic procedures.