

Anti-MSH2 monoclonal Antibody

Catalog Number: M00140-4


Overview

Product Name	Anti-MSH2 monoclonal Antibody
Reactive Species	Human
Description	Boster Bio Anti-MSH2 monoclonal Antibody catalog # M00140-4. Tested in WB,ICC applications. This antibody reacts with Human.
Application	ICC, WB
Clonality	Monoclonal 5D11
Formulation	1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Mouse
Uniprot ID	P43246


Technical Details

Immunogen	Recombinant full length Human MSH2.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1000 ICC: 1:50-200

Anti-MSH2 monoclonal Antibody (M00140-4) Images



Western blot detection of MSH2 antibody in HeLa, A431, A549 and 293T cell lysates using MSH2 antibody (1:500 diluted) .



Immunocytochemistry staining of HeLa cells using anti-MSH2 antibody (dilution 1:100).

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-MSH2 monoclonal Antibody

For Research Use Only. Not for use in diagnostic procedures.