

Anti-Cyclin D1 Rabbit Monoclonal Antibody

Catalog Number: M00149-5


Overview

Product Name	Anti-Cyclin D1 Rabbit Monoclonal Antibody
Reactive Species	Human, Mouse, Rat
Description	This Recombinant Rabbit Monoclonal Antibody detects endogenous levels of Cyclin D1 protein. Validated for research with WB,IHC,IF,IP,ELISA applications. Anti Cyclin D1 is reactive for Human, Mouse, Rat samples
Application	ELISA, IP, IF, IHC, WB
Clonality	Monoclonal
Formulation	PBS, 50% glycerol, 0.05% Proclin 300, 0.05% stabilizing protein This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C/1 year
Host	Rabbit
Uniprot ID	P24385


Technical Details

Isotype	IgG, Kappa
Purification	Protein A
Suggested Dilutions	IHC 1:200-1:1000 WB 1:1000-1:5000 IF 1:200-1:1000 ELISA 1:5000-1:20000 IP 1:50-1:200


Anti-Cyclin D1 Rabbit Monoclonal Antibody (M00149-5) Images




Human esophagus was stained with anti-Cyclin D1 rabbit antibody




Human tonsil was stained with anti-Cyclin D1 rabbit antibody




Rat spleen was stained with anti-Cyclin D1 rabbit antibody



Mouse colon was stained with anti-Cyclin D1 rabbit antibody



Mouse spleen was stained with anti-Cyclin D1 rabbit antibody



Various whole cell lysates were separated by 4-20% SDS-PAGE, and the membrane was blotted with anti-Cyclin D1 antibody. The HRP-conjugated Goat anti-Rabbit IgG(H + L) antibody was used to detect the antibody. Lane 1: Mouse spleen Lane 2: Rat spleen Lane 3: LnCap Lane 4: MCF7
Predicted band size: 34kDa Observed band size: 36kDa

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Cyclin D1 Rabbit Monoclonal Antibody

For Research Use Only. Not for use in diagnostic procedures.