

Anti-CD19 Reference Antibody (obexelimab)

Catalog Number: M00154-12


Overview

Product Name	Anti-CD19 Reference Antibody (obexelimab)
Reactive Species	Cynomolgus monkey, Human
Description	Boster Bio Anti-CD19 Reference Antibody (obexelimab) (Catalog # M00154-12). Tested in Flow Cytometry, ELISA, FTA. This antibody reacts with Cynomolgus monkey, Human. Endotoxin: < 0.227EU/ug,determined by LAL method. Expression system: CHO Cell
Application	ELISA, Flow Cytometry, Functional Assay, Kinetics
Clonality	Monoclonal
Formulation	Purified monoclonal antibody supplied in PBS, pH6.0, without preservative.This antibody is purified through a protein A column.
Storage Instructions	Maintain refrigerated at 2-8°C for up to 2 weeks. For long-term storage, aliquot and store at -20°C to avoid repeated freeze-thaw cycles.
Uniprot ID	P15391


Technical Details

Isotype	IgG1
Form	Liquid

Anti-CD19 Reference Antibody (obexelimab) (M00154-12) Images



Anti-CD19 Reference Antibody (obexelimab) on SDS-PAGE under reducing (R) condition. The gel was stained with Coomassie Blue. The purity of the protein is greater than 90%



The purity of Anti-CD19 Reference Antibody (obexelimab) is more than 95%

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-CD19 Reference Antibody (obexelimab)

For Research Use Only. Not for use in diagnostic procedures.