

Anti-ICAM1 / CD54 Reference Antibody (bersanlimab)

Catalog Number: M00171-7


Overview

Product Name	Anti-ICAM1 / CD54 Reference Antibody (bersanlimab)
Reactive Species	Human
Description	Boster Bio Anti-ICAM1 / CD54 Reference Antibody (bersanlimab) (Catalog # M00171-7). Tested in Flow Cytometry, ELISA, FTA. This antibody reacts with Human. Endotoxin: < 0.243EU/ug,determined by LAL method. Expression system: CHO Cell
Application	ELISA, Flow Cytometry, Functional Assay, Kinetics
Clonality	Monoclonal
Formulation	Purified monoclonal antibody supplied in PBS, pH6.0, without preservative.This antibody is purified through a protein A column.
Storage Instructions	Maintain refrigerated at 2-8°C for up to 2 weeks. For long-term storage, aliquot and store at -20°C to avoid repeated freeze-thaw cycles.
Uniprot ID	P05362


Technical Details

Isotype	IgG1
Form	Liquid

Anti-ICAM1 / CD54 Reference Antibody (bersanlimab) (M00171-7) Images



Anti-ICAM1/CD54 Reference Antibody (bersanlimab) on SDS-PAGE under reducing (R) condition. The gel was stained with Coomassie Blue. The purity of the protein is greater than 95.0%



The purity of Anti-ICAM1/CD54 Reference Antibody (bersanlimab) is more than 100%

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-ICAM1 / CD54 Reference Antibody (bersanlimab)

For Research Use Only. Not for use in diagnostic procedures.