

## Anti-PD1 monoclonal PDCD1 Antibody

Catalog Number: M00178-5


### Overview

Product Name	Anti-PD1 monoclonal PDCD1 Antibody
Reactive Species	Human
Description	Boster Bio Anti-PD1 monoclonal PDCD1 Antibody catalog # M00178-5. Tested in WB,ICC,IHC,Flow Cytometry applications. This antibody reacts with Human.
Application	Flow Cytometry, IHC, ICC, WB
Clonality	Monoclonal 3B12
Formulation	Mouse IgG1, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Mouse
Uniprot ID	Q15116


### Technical Details

Immunogen	Recombinant protein
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	ProA affinity purified
Suggested Dilutions	WB: 1:500-1:2,000 ICC: 1:50-1:200 IHC: 1:100-1:500 FC: 1:50-1:100

## Anti-PD1 monoclonal PDCD1 Antibody (M00178-5) Images



Western blot analysis of PD1 on PD1 transfected 293FT cell lysate using anti-PD1 antibody at 1/1,000 dilution.



Immunohistochemical analysis of paraffin-embedded human tonsil tissue using anti-PD1 antibody. Counter stained with hematoxylin.

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-PD1 monoclonal PDCD1 Antibody

For Research Use Only. Not for use in diagnostic procedures.