

## Anti-PDCD1 / PD-1 / CD279 Reference Antibody (nivolumab)

Catalog Number: M00178-8

### Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.




### Overview

|                      |                                                                                                                                                                                                                                                |
|----------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Name         | Anti-PDCD1 / PD-1 / CD279 Reference Antibody (nivolumab)                                                                                                                                                                                       |
| Reactive Species     | Human                                                                                                                                                                                                                                          |
| Description          | Boster Bio Anti-PDCD1 / PD-1 / CD279 Reference Antibody (nivolumab) (Catalog # M00178-8). Tested in Flow Cytometry, ELISA, FTA. This antibody reacts with Human. Endotoxin: < 0.245EU/ug,determined by LAL method. Expression system: CHO Cell |
| Application          | ELISA, Flow Cytometry, Functional Assay, Kinetics                                                                                                                                                                                              |
| Clonality            | Monoclonal                                                                                                                                                                                                                                     |
| Formulation          | Purified monoclonal antibody supplied in PBS, pH6.0, without preservative.This antibody is purified through a protein A column.                                                                                                                |
| Storage Instructions | Maintain refrigerated at 2-8°C for up to 2 weeks. For long-term storage, aliquot and store at -20°C to avoid repeated freeze-thaw cycles.                                                                                                      |
| Uniprot ID           | Q15116                                                                                                                                                                                                                                         |


### Technical Details

|         |        |
|---------|--------|
| Isotype | IgG4   |
| Form    | Liquid |


## Anti-PDCD1 / PD-1 / CD279 Reference Antibody (nivolumab) (M00178-8) Images



Anti-PDCD1/PD-1/CD279 Reference Antibody (nivolumab) on SDS-PAGE under reducing (R) condition. The gel was stained with Coomassie Blue. The purity of the protein is greater than 95%



The purity of Anti-PDCD1/PD-1/CD279 Reference Antibody (nivolumab) is more than 99.03%



Immobilized human PD 1 His at 2 µg/mL can bind Anti-PDCD1/PD-1/CD279 Reference Antibody (nivolumab)

[Anti-PDCD1 / PD-1 / CD279 Reference Antibody \(nivolumab\)](#)

For Research Use Only. Not for use in diagnostic procedures.