

Anti-Caveolin-1 CAV1 Rabbit Monoclonal Antibody

Catalog Number: M00179


Overview

Product Name	Anti-Caveolin-1 CAV1 Rabbit Monoclonal Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-Caveolin-1 CAV1 Rabbit Monoclonal Antibody catalog # M00179. Tested in WB, IHC, ICC/IF, Flow Cytometry applications. This antibody reacts with Human, Mouse, Rat.
Application	Flow Cytometry, IF, IHC, ICC, WB
Clonality	Monoclonal GE-3
Formulation	Rabbit IgG in stabilizing components, phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. *This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q03135


Technical Details

Immunogen	A synthesized peptide derived from human Caveolin-1
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5mg/ml
Purification	Affinity-chromatography
Suggested Dilutions	WB 1:500-2000 IHC 1:50-200 ICC/IF 1:50-200 FC 1:20


Anti-Caveolin-1 CAV1 Rabbit Monoclonal Antibody (M00179) Images




All lanes use the Antibody at 1:2W dilution for 1 hour at room temperature.



Immunofluorescent analysis of A431 cells, using Caveolin-1 Antibody .




Western blot analysis of Caveolin-1 expression in A431 cell lysate.



All lanes use the Antibody at 1:2W dilution for 1 hour at room temperature.

Immunohistochemical analysis of paraffin-embedded mouse lung, using Caveolin-1 Antibody.



1 Publications Citing This Product

1. PubMed ID: 24830555, MicroRNA-199a-5p Affects Porcine Preadipocyte Proliferation and Differentiation

Visit bosterbio.com/anti-caveolin-1-rabbit-monoclonal-antibody-m00179-boster.html to see all 1 publications.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Caveolin-1 CAV1 Rabbit Monoclonal Antibody

For Research Use Only. Not for use in diagnostic procedures.