

Anti-ITCH Rabbit Monoclonal Antibody

Catalog Number: M00195-2

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Overview

Product Name	Anti-ITCH Rabbit Monoclonal Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-ITCH Rabbit Monoclonal Antibody catalog # M00195-2. Tested in WB, IHC applications. This antibody reacts with Human, Mouse, Rat.
Application	IHC, WB
Clonality	Monoclonal 23I97
Formulation	Rabbit IgG in stabilizing components, phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. *This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q96J02


Technical Details

Immunogen	A synthesized peptide derived from human ITCH
Isotype	IgG
Form	Liquid
Concentration	0.5mg/ml
Purification	Affinity-chromatography

Suggested Dilutions

WB 1:1000-5000
IHC 1:50-200

Anti-ITCH Rabbit Monoclonal Antibody (M00195-2) Images



Western blot analysis of ITCH/AIP4 expression in K562 cell lysate.

[Anti-ITCH Rabbit Monoclonal Antibody](#)

For Research Use Only. Not for use in diagnostic procedures.