

Anti-CD80/B7 1 Rabbit Monoclonal Antibody

Catalog Number: M00196


Overview

| | |
|----------------------|--|
| Product Name | Anti-CD80/B7 1 Rabbit Monoclonal Antibody |
| Reactive Species | Human |
| Description | Boster Bio Anti-CD80/B7 1 Rabbit Monoclonal Antibody catalog # M00196. Tested in WB, IHC applications. This antibody reacts with Human. |
| Application | IHC, WB |
| Clonality | Monoclonal AOEH-3 |
| Formulation | Rabbit IgG in stabilizing components, phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. *This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required. |
| Storage Instructions | Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles. |
| Host | Rabbit |
| Uniprot ID | P33681 |


Technical Details

| | |
|---------------------|---|
| Immunogen | A synthesized peptide derived from human CD80 |
| Isotype | Rabbit IgG |
| Form | Liquid |
| Concentration | 0.5mg/ml |
| Purification | Affinity-chromatography |
| Suggested Dilutions | WB 1:500-2000 IHC 1:50-200 |

Anti-CD80/B7 1 Rabbit Monoclonal Antibody (M00196) Images



All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.



Western blot analysis of CD80 expression in Raji cell lysate.

1 Publications Citing This Product

1. PubMed ID: 33537094, Chen X,Gao A,Zhang F,Yang Z,Wang S,Fang Y,Li J,Wang J,Shi W,Wang L, Zheng Y,Sun Y.ILT4 inhibition prevents TAM- and dysfunctional T cell-mediated immunosuppression and enhances the efficacy of anti-PD-L1 therapy in NSCLC with EGFR activation.Theranostics.

Visit bosterbio.com/anti-cd80-rabbit-monoclonal-antibody-m00196-boster.html to see all 1 publications.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-CD80/B7 1 Rabbit Monoclonal Antibody

For Research Use Only. Not for use in diagnostic procedures.