

Anti-FGFR3 Rabbit Monoclonal Antibody

Catalog Number: M00200-1


Overview

Product Name	Anti-FGFR3 Rabbit Monoclonal Antibody
Reactive Species	Human
Description	Boster Bio Anti-FGFR3 Rabbit Monoclonal Antibody catalog # M00200-1. Tested in WB, IHC, ICC/IF, Flow Cytometry applications. This antibody reacts with Human.
Application	Flow Cytometry, IF, IHC, ICC, WB
Clonality	Monoclonal AACG-6
Formulation	Rabbit IgG in stabilizing components, phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. *This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P22607


Technical Details

Immunogen	A synthesized peptide derived from human FGFR3
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5mg/ml
Purification	Affinity-chromatography
Suggested Dilutions	WB 1:500-2000 IHC 1:20-100 ICC/IF 1:50~1:200 FC 1:50


Anti-FGFR3 Rabbit Monoclonal Antibody (M00200-1) Images



Western blot analysis of FGFR3 using anti-FGFR3 antibody (M00200-1). Electrophoresis was performed on a 5-20% SDS-PAGE gel at 70V (Stacking gel) / 90V (Resolving gel) for 2-3 hours. The sample well of each lane was loaded with 30 ug of sample under reducing conditions. Lane 1: human Hela whole cell lysates, Lane 2: human HepG2 whole cell lysates, Lane 3: human RT4 whole cell lysates, Lane 4: human CACO-2 whole cell lysates, Lane 5: human HEL whole cell lysates, Lane 6: human 293T whole cell lysates. After electrophoresis, proteins were transferred to a nitrocellulose membrane at 150 mA for 50-90 minutes. Blocked the membrane with 5% non-fat milk/TBS for 1.5 hour at RT. The membrane was incubated with rabbit anti-FGFR3 antigen affinity purified monoclonal antibody (M00200-1) at 1:500 overnight at 4°C, then washed with TBS-0.1%Tween 3 times with 5 minutes each and probed with a goat anti-rabbit IgG-HRP secondary antibody at a dilution of 1:5000 for 1.5 hour at RT. The signal is developed using an Enhanced Chemiluminescent detection (ECL) kit (Catalog # EK1002) with Tanon 5200 system. A specific band was detected for FGFR3 at approximately 125 kDa. The expected band size for FGFR3 is at 100 kDa.




All lanes use the Antibody at 1:3K dilution for 1 hour at room temperature.



IHC analysis of FGFR3 using anti-FGFR3 antibody (M00200-1). FGFR3 was detected in a paraffin-embedded section of human bladder cancer tissue. Heat mediated antigen retrieval was performed in EDTA buffer (pH 8.0, epitope retrieval solution). The tissue section was blocked with 10% goat serum. The tissue section was then incubated with 1:50 rabbit anti-FGFR3 Antibody (M00200-1) overnight at 4°C. Peroxidase Conjugated Goat Anti-rabbit IgG was used as secondary antibody and incubated for 30 minutes at 37°C. The tissue section was developed using HRP Conjugated Rabbit IgG Super Vision Assay Kit (Catalog # SV0002) with DAB as the chromogen.

IHC analysis of FGFR3 using anti-FGFR3 antibody (M00200-1). FGFR3 was detected in a paraffin-embedded section of human breast cancer tissue. Heat mediated antigen retrieval was performed in EDTA buffer (pH 8.0, epitope retrieval solution). The tissue section was blocked



with 10% goat serum. The tissue section was then incubated with 1:50 rabbit anti-FGFR3 Antibody (M00200-1) overnight at 4°C. Peroxidase Conjugated Goat Anti-rabbit IgG was used as secondary antibody and incubated for 30 minutes at 37°C. The tissue section was developed using HRP Conjugated Rabbit IgG Super Vision Assay Kit (Catalog # SV0002) with DAB as the chromogen.

1 Publications Citing This Product

1. PubMed ID: 33626078, Nannapaneni S,Griffith CC,Magliocca KR,Chen W,Lyu X,Chen Z,Wang D,Wang X,Shin DM,Chen ZG,Saba NF.Co-expression of fibroblast growth factor receptor 3 with mutant p53, and its association with worse outcome in oropharyngeal squamous cell carcinoma.PLoS One

Visit bosterbio.com/anti-fgfr3-rabbit-monoclonal-antibody-m00200-1-boster.html to see all 1 publications.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-FGFR3 Rabbit Monoclonal Antibody

For Research Use Only. Not for use in diagnostic procedures.