

## Anti-Parathyroid hormone PTH Monoclonal Antibody

Catalog Number: M00208-1


### Overview

Product Name	Anti-Parathyroid hormone PTH Monoclonal Antibody
Reactive Species	Human
Description	Boster Bio Anti-Parathyroid hormone PTH Monoclonal Antibody catalog # M00208-1. Tested in ELISA, WB applications. This antibody reacts with Human.
Application	ELISA, WB
Clonality	Monoclonal 2C7
Formulation	Ascitic fluid containing 0.03% sodium azide.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Mouse
Uniprot ID	P01270

### Technical Details

Immunogen	Purified recombinant fragment of human PTH (aa1-115) expressed in E. Coli.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	Affinity purification
Suggested Dilutions	WB 1:500-1:2000 ELISA 1:10000

## Anti-Parathyroid hormone PTH Monoclonal Antibody (M00208-1) Images



Western blotting validation for Anti-Parathyroid hormone PTH Monoclonal Antibody M00208-1 Western Blot (WB) analysis using PTH Monoclonal Antibody against PTH-hlgGfC transfected HEK293 cell lysate. Electrophoresis was performed on a SDS-PAGE gel. To determine SDS-PAGE gel concentration

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Parathyroid hormone PTH Monoclonal Antibody

For Research Use Only. Not for use in diagnostic procedures.