

## Anti-Parathyroid Hormone Rabbit Monoclonal Antibody

Catalog Number: M00208-2


### Overview

|                      |   |
|----------------------|---|
| Product Name         | Anti-Parathyroid Hormone Rabbit Monoclonal Antibody   |
| Reactive Species     | Human   |
| Description          | Boster Bio Anti-Parathyroid Hormone Rabbit Monoclonal Antibody catalog # M00208-2. Tested in WB, IHC applications. This antibody reacts with Human. |
| Application          | IHC, WB   |
| Clonality            | Monoclonal 17P72  |
| Formulation          | Rabbit IgG in phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol, 0.4-0.5mg/ml BSA.                                 |
| Storage Instructions | Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.          |
| Host                 | Rabbit  |
| Uniprot ID           | P01270  |

### Technical Details

|                     |   |
|---------------------|---|
| Immunogen           | A synthesized peptide derived from human Parathyroid Hormone  |
| Isotype             | IgG   |
| Form                | Liquid  |
| Concentration       | Actual concentration vary by lot. Use suggested dilution ratio to decide dilution procedure.  |
| Purification        | Affinity-chromatography   |
| Suggested Dilutions | <p>Dilute the sample so that the expected range of concentrations fall within the detection range of this kit.</p> <p>If the expected range of concentration is unknown, a pilot test should be conducted to decide the optimal dilution ratio for your samples.</p> <p>Some PubMed article(s) citing the expression level of this target are as follows:<br/>Boster Bio's internal QC testing used WB 1:500-1:2000<br/>IHC 1:50-1:200</p> <p>For protocols, please visit <a href="https://www.bosterbio.com/protocol-and-troubleshooting/">https://www.bosterbio.com/protocol-and-troubleshooting/</a></p> |

## Anti-Parathyroid Hormone Rabbit Monoclonal Antibody (M00208-2) Images



Western blot analysis of Parathyroid Hormone expression in human Parathyroid Hormone recombinant protein.

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Parathyroid Hormone Rabbit Monoclonal Antibody