

Anti-SNCA Antibody

Catalog Number: M00215-4

About SNCA

May be involved in the regulation of dopamine release and transport. Induces fibrillization of microtubule-associated protein tau. Reduces neuronal responsiveness to various apoptotic stimuli, leading to a decreased caspase-3 activation.

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Overview


Product Name	Anti-SNCA Antibody
Reactive Species	Human, Mouse
Description	Boster Bio Anti-SNCA Antibody (Catalog # M00215-4). Tested in WB application(s). This antibody reacts with Human, Mouse.
Application	WB
Clonality	Monoclonal 1853CT506.24.16
Formulation	Purified monoclonal antibody supplied in PBS with 0.09% (W/V) sodium azide.
Storage Instructions	Maintain refrigerated at 2-8°C for up to 2 weeks. For long-term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Host	Mouse
Uniprot ID	P37840

Technical Details

Immunogen	This antibody is generated from a mouse immunized a recombinant protein from human.
Predicted Reactive Species	Human, Mouse, Rat
Isotype	IgG1,kappa
Purification	This antibody is purified through a protein G column, followed by dialysis against PBS.
Suggested Dilutions	WB: 1:2000 Tested Species: In-house tested species with positive results. Other applications have not been tested. Optimal dilutions

	should be determined by end users.
--	------------------------------------

Anti-SNCA Antibody (M00215-4) Images



All lanes : Anti-SNCA at 1:2000 dilution
Lane 1: Human brain lysate
Lane 2: Mouse brain lysate
Lysates/proteins at 20 μ g per lane. Secondary Goat Anti-mouse IgG, (H+L), Peroxidase conjugated at 1/10000 dilution. Predicted band size : 14 kDa
Blocking/Dilution buffer: 5% NFDM/TBST.

Anti-SNCA Antibody

For Research Use Only. Not for use in diagnostic procedures.