

## Anti-LKB1 STK11 Monoclonal Antibody

Catalog Number: M00217

### Overview

Product Name	Anti-LKB1 STK11 Monoclonal Antibody
Reactive Species	Human, Monkey, Mouse
Description	Boster Bio Anti-LKB1 STK11 Monoclonal Antibody catalog # M00217. Tested in ELISA, Flow Cytometry, WB applications. This antibody reacts with Human, Monkey, Mouse.
Application	ELISA, Flow Cytometry, WB
Clonality	Monoclonal 4H12
Formulation	Ascitic fluid containing 0.03% sodium azide.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Mouse
Uniprot ID	Q15831

### Technical Details

Immunogen	Purified recombinant fragment of human LKB1 expressed in E. Coli.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	Affinity purification
Suggested Dilutions	WB 1:500-1:2000 FC 1:200-1:400 ELISA 1:10000

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-LKB1 STK11 Monoclonal Antibody

For Research Use Only. Not for use in diagnostic procedures.