

Anti-AXL / UFO Reference Antibody (enapotamab vedotin)

Catalog Number: M00226-4


Overview

Product Name	Anti-AXL / UFO Reference Antibody (enapotamab vedotin)
Reactive Species	Cynomolgus monkey, Human
Description	Boster Bio Anti-AXL / UFO Reference Antibody (enapotamab vedotin) (Catalog # M00226-4). Tested in Flow Cytometry, ELISA, FTA. This antibody reacts with Cynomolgus monkey, Human. Endotoxin: < 0.298EU/ug,determined by LAL method. Expression system: CHO Cell
Application	ELISA, Flow Cytometry, Functional Assay, Kinetics
Clonality	Monoclonal
Formulation	Purified monoclonal antibody supplied in PBS, pH6.0, without preservative.This antibody is purified through a protein A column.
Storage Instructions	Maintain refrigerated at 2-8°C for up to 2 weeks. For long-term storage, aliquot and store at -20°C to avoid repeated freeze-thaw cycles.
Uniprot ID	P30530

Technical Details

Isotype	IgG1
Form	Liquid

Anti-AXL / UFO Reference Antibody (enapotamab vedotin) (M00226-4) Images



The purity of Anti-AXL/UFO Reference Antibody (enapotamab vedotin) is more than 99.9%

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-AXL / UFO Reference Antibody (enapotamab vedotin)

For Research Use Only. Not for use in diagnostic procedures.