

Anti-GAPDH Rabbit Monoclonal Antibody Biotin Conjugated

Catalog Number: M00227-Biotin

Overview

Product Name	Anti-GAPDH Rabbit Monoclonal Antibody Biotin Conjugated
Reactive Species	Human, Mouse, Rat
Application	WB, IHC, ELISA
Clonality	Monoclonal BH-7
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% Na ₃ N.
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	P04406

Technical Details

Immunogen	A synthesized peptide derived from human GAPDH
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Affinity-chromatography
Conjugate	Biotin
Suggested Dilutions	Western blot, Optimal dilutions should be determined by end users. Immunohistochemistry (Paraffin-embedded Section), Optimal dilutions should be determined by end users. ELISA, Optimal dilutions should be determined by end users.

153 Publications Citing This Product

1. PubMed ID: 10.3892/etm.2018.5822, Effective robotic assistive pattern of treadmill training for spinal cord injury in a rat model
2. PubMed ID: 31173299, Wang J,Fang Y,Liu YF,Wang X,Wang XL,Wang RY,Meng ZD.MiR-154 inhibits cells proliferation and metastasis in melanoma by targeting AURKA and serves as a novel prognostic indicator.Eur Rev Med Pharmacol Sci.2019 May;23(10):4275-4284.doi:10.26355/eurrev_201905_17932.PMID:31173299.

3. PubMed ID: 31058194, Duan Y, Tan Z, Yang M, Li J, Liu C, Wang C, Zhang F, Jin Y, Wang Y, Zhu L. PC-3-Derived Exosomes Inhibit Osteoclast Differentiation by Downregulating miR-214 and Blocking NF-kappaB Signaling Pathway. Biomed Res Int. 2019 Apr 1; 2019:8650846. doi:10.1155/2019/8650846. PMID:31058194; PMCID:PMC6463683.

Visit bosterbio.com/anti-gapdh-rabbit-monoclonal-antibody-m00227-boster.html to see all 153 publications.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-GAPDH Rabbit Monoclonal Antibody - Biotin

For Research Use Only. Not for use in diagnostic procedures.