

## Anti-CD147 Rabbit Monoclonal Antibody

Catalog Number: M00248-7


### Overview

Product Name	Anti-CD147 Rabbit Monoclonal Antibody
Reactive Species	Mouse
Description	Boster Bio Anti-CD147 Rabbit Monoclonal Antibody catalog # M00248-7. Tested in WB, IHC, ICC/IF, IP applications. This antibody reacts with Mouse.
Application	IP, IF, IHC, ICC, WB
Clonality	Monoclonal 30B80
Formulation	Rabbit IgG in stabilizing components, phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. *This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P18572

### Technical Details

Immunogen	A synthesized peptide derived from human CD147
Isotype	IgG
Form	Liquid
Concentration	0.5mg/ml
Purification	Affinity-chromatography
Suggested Dilutions	WB 1:500-1:2000 IHC 1:50-1:200 ICC/IF 1:50-1:200 IP 1:50

## Anti-CD147 Rabbit Monoclonal Antibody (M00248-7) Images



Western blot analysis of CD147 expression in RAW 264.7 cell lysate.

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-CD147 Rabbit Monoclonal Antibody

For Research Use Only. Not for use in diagnostic procedures.