

## Anti-Prealbumin TTR Rabbit Monoclonal Antibody

Catalog Number: M00290


### Overview

|                      |                                                                                                                                                                                                                                                                                                                                                                                                                                                    |
|----------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Name         | Anti-Prealbumin TTR Rabbit Monoclonal Antibody                                                                                                                                                                                                                                                                                                                                                                                                     |
| Reactive Species     | Human                                                                                                                                                                                                                                                                                                                                                                                                                                              |
| Description          | Boster Bio Anti-Prealbumin TTR Rabbit Monoclonal Antibody catalog # M00290. Tested in WB, IHC, IP applications. This antibody reacts with Human.                                                                                                                                                                                                                                                                                                   |
| Application          | IP, IHC, WB                                                                                                                                                                                                                                                                                                                                                                                                                                        |
| Clonality            | Monoclonal AOFI-20                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| Formulation          | Rabbit IgG in stabilizing components, phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.<br>*This antibody is supplied in a stabilized formulation.<br>Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required. |
| Storage Instructions | Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.                                                                                                                                                                                                                                                                                                         |
| Host                 | Rabbit                                                                                                                                                                                                                                                                                                                                                                                                                                             |
| Uniprot ID           | P02766                                                                                                                                                                                                                                                                                                                                                                                                                                             |

### Technical Details

|                     |                                                     |
|---------------------|-----------------------------------------------------|
| Immunogen           | A synthesized peptide derived from human Prealbumin |
| Isotype             | Rabbit IgG                                          |
| Form                | Liquid                                              |
| Concentration       | 0.5mg/ml                                            |
| Purification        | Affinity-chromatography                             |
| Suggested Dilutions | WB 1:500-2000<br>IHC 1:50-200<br>IP 1:20            |

## Anti-Prealbumin TTR Rabbit Monoclonal Antibody (M00290) Images



Western blot analysis of Prealbumin gamma expression in human fetal kidney lysate.

### 1 Publications Citing This Product

1. PubMed ID: 30340515, Expression of renal cell markers and detection of 3p loss links endolymphatic sac tumor to renal cell carcinoma and warrants careful evaluation to avoid %u2026

Visit [bosterbio.com/anti-prealbumin-rabbit-monoclonal-antibody-m00290-boster.html](http://bosterbio.com/anti-prealbumin-rabbit-monoclonal-antibody-m00290-boster.html) to see all 1 publications.

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Prealbumin TTR Rabbit Monoclonal Antibody

For Research Use Only. Not for use in diagnostic procedures.